Replens Long Lasting Vaginal Moisturiser x 12 Pre-filled Applicators



Category

Female Health and Genitourinary Manufacturer **CHURCH & DWIGHT UK** LTD Pack size 12 RRP £11.54 Legal Status MD

0 Uses

Long lasting vaginal moisturiser. Helps replenish internal vaginal moisture, alleviate vaginal dryness symptoms such as discomfort, itch, irritation and provide more comfortable intimacy. One application provides up to 3 days' relief.



Format and Ingredients

Pre-filled applicator

Product data © Communications International Group

Dosage

Adults: apply once every 3 days. Remove the applicator from the sealed wrapper and grip firmly by the thick end. Shake down like a thermometer to ensure that the contents are at the thin end. Bend the thin pointed end of the applicator back and forth repeatedly. Snap off the thin pointed end of the applicato and discard the break-away tip. The applicator may be inserted into the vagina while you are in a sitting position or when lying on your back with knees bent. Gently insert the thin end well into the vagina as deeply as it will go comfortably. Squeeze the thick end of the applicator firmly to deposit gel.